8
1L

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PHILIPS

Technical data

Power supply		Optimus 50	Optimus 65	Optimus 80
Mains voltage	400 V (±10%); 50 Hz or 60 Hz, 3-phase If the generators are connected up to mains which is not 400 V (±10%), a mains adaptation transformer is required.			
Max. mains resistance at 400 V	0,3 Ω	0,2 Ω	0,2 Ω	
Max. mains current at 400 V	145 A	190 A	230 A	
Fuses	35 A (slow)	50 A (slow)	50 A (slow)	
X-ray generator				
High-voltage generator	The converter generator generates high voltage equivalent to DC voltage			
Nominal power (IEC)	50 kW	65 kW	80 kW	
Max. tube voltage	150 kV	150 kV	150 kV	
Max. tube current (at 70 kV)	650 mA	900 mA	1100 mA	
Max. frame rates	6 fps for camera technique and DSI			
Compatible X-ray tubes	Up to 3 double-focus tubes with a total of 6 different focal spots can be used. SRO- and RO tube assemblies are compatible.			
Tube load protection	The tube is thermally monitored by the tube load protection system. If the load is high, tube output is reduced; exposures being made are not aborted. The tube load display on the generator control desk indicates the tube power available.			
Nominal tube output				
Ratings	650 mA at 70 kV	900 mA at 70 kV	1100 mA at 70 kV	
	625 mA at 80 kV	812 mA at 80 kV	1000 mA at 80 kV	
	500 mA at 100 kV	650 mA at 100 kV	800 mA at 100 kV	
	400 mA at 125 kV	520 mA at 125 kV	640 mA at 125 kV	
	333 mA at 150 kV	433 mA at 150 kV	533 mA at 150 kV	
Ambient conditions				
Temperature	+10 °C to +40 °C			
Rel. humidity	15% to 90%, without condensation			
Rel. atmospheric pressure	50 kPa to 110 kPa			
Sound emission	55 dBa max., averaged over 1 h, at a distance of 1 m			
Ventilation demand	2.5 l/s max.			

Radiography and Fluoroscopy

Radiography	Optimus 50	Optimus 65	Optimus 80
Exposure techniques:			
One-factor technique	kV, continuously falling load, AMPLIMAT ¹⁾		
Two-factor technique	kV, mA, constant load, AMPLIMAT ¹⁾		
Two-factor technique	kV, mAs, constant load		
Three-factor technique	kV, mA, s, constant load		
Three-factor technique	kV, mAs, s, constant load		

¹⁾ AMPLIMAT-Automatic exposure control

- up to 5 ionization measuring chambers or
- 4 ionization measuring chambers and 1 photomultiplier (secondary electron multiplier) or
- 3 ionization measuring chambers, 1 photomultiplier and 1 photodiode

Setting ranges

Tube voltage	40 kV to 150 kV programmable in: 27 large steps or 111 small steps of 1 kV each		
Radiography without AMPLIMAT:			
Tube current	1 mA to 650 mA ²⁾	1 mA to 900 mA ²⁾	1 mA to 1100 mA ²⁾
mAs product	0.5 mAs to 850 mAs		
Steps	For tube current and mAs product, programmable as required in: 6.25%, 12.5% or 25% steps		
Exposure times	<ul style="list-style-type: none">• kV-mA-s-technique: you select in steps from 1 ms to 6 s (16 s ³⁾)• kV-mAs-s-technique: from 1 ms to 6 s the generator switches automatically		
Radiography with AMPLIMAT:			
mAs product	0.5 mAs to 800 mAs ²⁾ , the maximum can be limited to 600 mAs		
Switch-off dose	0.5 µGy to 100 µGy		
Exposure	1 ms to 4 s or to 6 s with TDC (Tomographic Density Control) The display panel on the control desk shows the switched mAs- and ms values after exposure.		
Screen-film combinations	From a library of approx. 100 film and screen data, 5 screen-film combinations can be read in. The values of new combinations can be integrated by Customer Service at any time.		

²⁾ The ultimate value depends on the particular tube nomogram.

³⁾ For special application techniques.

Fluoroscopy

Image intensifier TV system	1
Tube voltage	40 kV to 110 kV (programmable up to 125 kV on request)
Tube current	0.1 mA to 6 mA (maximum selectable between 3 mA and 6 mA)

Quick fluoroscopy

Tube voltage	40 kV to 110 kV (programmable up to 125 kV on request)
Tube current	0.1 mA to 6 mA (maximum selectable between 3 mA and 6 mA)

Adjustment times:


Standard fluoroscopy	0.5 s to 1.0 s (dependent on absorption)
Quick fluoroscopy	0.2 s to 0.4 s (dependent on absorption)

Tube description

The many years of experience at Philips in manufacturing Roralix X-ray tubes led to the use of state-of-the-art manufacturing technologies and materials. For example, the all-metal rotating anode has high permissible temperatures so heat dissipation is improved and the tube cools down quickly.

The Roralix X-ray tube consists of the glass vacuum envelope containing the following main components:

- cathode (K),
- rotating anode (A) and
- rotor with ball bearings (R).

Precision and quality

All Roralix X-ray tubes are subject of rigorous individual testing under practical conditions before they leave the factory. In this way Philips ensures that all the ratings are met and a long tube life is ensured. Radiation quality is assured and radiation protection is guaranteed.

Anode and cathode

Fitting the anode and cathode into the glass envelope calls for sophisticated production techniques and skilful craftsmanship. This precision work is performed by a team of specialists with years of experience. It is an area where the specialised know-how of glassblowers still counts.

Precision adjustment of a cathode

Technical data

	Roralix X-ray tubes		
	RO 12 30	RO 17 50	RO 3050 Rc
Nominal anode input power ¹⁾ [kW]	12	17	30
	30	50	50
Nominal focal spot values	0,6	0,6	1,2
	1,2	1,3	1,8
Anode angle	16°	13°	16°
Max. X-ray field at SID = 100 cm	53 x 53	43 x 43	53 x 53
Max. anode heat content			
— anode	200	220	200
— X-ray tube assembly	1260	1260	1260
Max. heat dissipation in continuous mode with rotating anode:			
— anode	350	450	350
— X-ray tube assembly:			
- with fan cooling	350	350	350
- without fan cooling	250	250	250
Nominal X-ray tube voltage	125 or 150	150	125 or 150
Cooling rate of anode	1300	1300	1300

¹⁾ based on 250 W equivalent anode input power